



Top-tier sterilization: performance and sustainability



#### Treating used water

The **recirculation filter** incorporated in the secondary reservoir allows used water to be demineralised and purified. This process makes it reusable in subsequent sterilization cycles and automatically controls the water delivery reservoir. Just 5 litres of mains water are needed to perform 50 sterilization cycles: a market-leading specific consumption of just 100 ml/cycle.



#### Water quality control

The built-in conductivity sensor ensures compliance with **water quality** parameters, protecting the internal components of the device. It also monitors the status of the filtration system. Clear messages on the display inform the user when the filters need replacing.



### Replacing the filters

Easy access to the reservoirs and carefully designed filter ergonomics make replacement easy.

After being informed by the on-screen messaging system that the filters need replacing, the user can follow the instructions contained in the preloaded video tutorials and proceed independently with replacement of integrated filters.

## Water savings

## 50 sterilization cycles with only 5 litres of mains water.

Uses normal mains water and can retrieve used water to employ it in subsequent sterilization cycles. These features minimise waste, reduce costs and significantly increase efficiency, eliminating any demineralised water supply-storage problems.





#### Integrated demineralisation

A demineralisation filter housed in the main tank lets users fill the reservoir with ordinary mains water, resulting in excellent rationalisation of the available space, streamlined processes and significant savings. Moreover, the filtration system delivers never-before-seen operational simplicity. To make water management even easier, optional automatic filling systems are available.



## Ease of use and high performance

The 7" touchscreen is as simple and user-friendly as a tablet.

Screen size and performance - plus meticulous design of the user interface - allow users to start sterilizing immediately with ease.

Video tutorials and luminous indicators put it in a class of its own.

The adopted components and the advanced control software provide outstanding performance.



#### LED bar

The device also communicates using light. The **LED bar incorporated** in the door provides immediate information on machine status by varying the colours according to the phase of operation. Consequently, users can see which phase of the process the autoclave is performing at a glance, even from some distance away.



# Shorter cycle times and Modular Drying

This sterilizer meets the needs of even the most demanding users. Reliable high-performance internal components and a control system that optimises each phase reduce sterilization cycle times. The exclusive Modular Drying function lets users select the load quantity to be sterilized. This adapts the duration of the drying phase, thus shortening cycle times. A complete B cycle in just 31 minutes.



## NFC - Near Field Communication

The **NFC sensor** simplifies the assistant's workflow. Activating user ID ensures only authorised personnel can access the machine functions. To activate, simply bring the NFC card, supplied with the sterilizer and associated with the user, up to the machine.



### Video tutorials, on-display info

High-resolution images make the 7" display a valuable tool for all users. Simple video tutorials provide important use and maintenance instructions. A full complement of information - from indications on load positioning to instructions for periodic maintenance - is readily accessible on the machine. What's more, cycle reports saved on the autoclave memory can be viewed on the screen.



# Customisable sterilization cycles

The autoclave has five sterilization cycles with pre-set parameters plus a **custom cycle** that can be personalised by the user. Additional functions allow the sterilizer to be used more effectively. The **delayed start** function lets you programme a deferred cycle start; the **extra drying** function is useful with particularly difficult loads as it lets you extend drying times.





#### Di.V.A.

**Di.V.A.** is an optional Cloud environment for sterilization devices. All saved data is made available to operators, allowing them to:

- · view reports of completed cycles
- $\cdot$  check device efficiency and status
- · access video tutorials and user manuals
- · consult machine usage statistics

# Connectivity and data management

Integrated Wi-Fi, Cloud platforms for remote assistance and a mobile APP for optimised management.

Connecting Supreme to the network via the integrated Wi-Fi system or Ethernet links the sterilizer to an integrated data management system, essential for advanced traceability. Users can also take advantage of other features such as automatically saved cycle reports, local or remote consultation and printer management.



#### **EasyCheck**

An optional in-Cloud platform that delivers remote technical support and significantly reduces service times, **EasyCheck** guarantees efficiency and reliability. Software updates, maintenance interval checks, remote verification of components and sensors are just some of the tools available to ensure long-lasting autoclave performance.



#### SterilConnect

A highly innovative tool users will appreciate. Any network-connected mobile device - tablet or smartphone - can be used to interact with the autoclave, simplifying daily management of the sterilization process. A clear yet complete graphic interface lets you monitor device status and manage the various sterilizer functions.



#### QR code labels

Activating this function lets users create adhesive QR code labels to apply to sterile packages at the end of each cycle. **The QR code** contains unique information that identifies the sterilizer, the cycle used, its outcome and the expiry date of the sterile package. Scanning the QR code lets you associate the data it contains with the patient, ensuring comprehensive traceability. It's also possible to print Barcode labels.



#### **Shared printer**

By producing labels or reports, the optional printer allows completion of the sterilization process, providing users with the basis for proper traceability. To optimise management, the **printer can be shared** by multiple sterilizers, where enabled, and connected to the same network. One fast, reliable printing station rationalises use of space and delivers significant savings.



## **Automated** sterilization protocol

The Mocom response to ensure correct implementation of sterilization protocols and protection of patients and health workers.









**TETHYS H10** 

resulting in unparalleled performance.

**THALYA - THALYA PLUS** 

The sterilization protocol begins with collecting instruments and materials.

undergo a cleaning, disinfection and drying process within a single cycle.

The Manual version delivers user-friendliness and outstanding reliability.





Disinfection Washing



Rinsing













**Packaging** 







**MILLSEAL** 



Designed to give you the best. The best performance in a simple and fast work-flow. An innovative thermal disinfector that replaces the numerous manual tasks typical of the stages preceding sterilization, thus reducing the workloads. Equipped with the optional HMD accessory, Tethys H10 extends the reconditioning process to rotating instruments,

Thalya is the dentist's ally for effective maintenance of rotary instruments. Simple and user-friendly, it allows perfect lubrication and purge maintenance of turbines and handpieces. In the Thalya Plus version, the instruments also

This thermal sealer range meets every possible need. The EVO version allows the automated creation of pouches. The Rolling version, also available with printer, combines sealing speed with flexible sizing of the packages to be sealed.





**MYTRACE** 







**STORAGE** 

If carried out correctly, storage ensures long-lasting the instrument sterility.















Technical data	Supreme 17 Supreme 22		Supreme 28			
Power supply	2220/240 V 50Hz 220/230 V 60Hz					
Rated power	2300 W					
External dimensions (L x H x D)	500 x 490 x 600 mm					
Chamber dimensions (Diam. x Depth)	250 x 350 mm	250 x 450 mm	280 x 450 mm			
Total weight	47 kg	50 kg	51 kg			
Tank capacity	61					
Autonomy	From 7 to 11 cycles	From 6 to 10 cycles	From 5 to 9 cycles			

Cycle type	Sterilization time	B Futura 17	B Futura 22	B Futura 28
,	(min.)	Cycle times inclu	ding sterilization tim time (min.)	nes, net of drying
В	4	39 (31)	42 (34)	52 (44)
В	20	53 (45)	59 (51)	63 (55)
S	4	21	24	27
S	4	32 (25)	37 (30)	41 (34)
В	18	53 (45)	57 (49)	67 (59)
Users can personalise the cycle with temperatures of 134 °C/121 °C, sterilization times starting from 4' (134 °C) or 20' (121 °C) and drying times from 5 to 30 min.				
	TEST	18	18	19
	TEST	20	24	28
	TEST	42	46	51
	type  B B S S	type time  (min.)  B 4 B 20 S 4 S 4 B 18 S Users can personal times starting from TEST TEST	type         time         B Futura 17           (min.)         Cycle times included and the properties of the cycle times included and the properties of	type         time         B Futura 17         B Futura 22           (min.)         Cycle times including sterilization tim time (min.)           B         4         39 (31)         42 (34)           B         20         53 (45)         59 (51)           S         4         21         24           S         4         32 (25)         37 (30)           B         18         53 (45)         57 (49)           S         Users can personalise the cycle with temperatures of 134 °C times starting from 4' (134 °C) or 20' (121 °C) and drying times           TEST         18         18           TEST         20         24

<sup>\*</sup> Drying varies according to machine model and volume Note: times in brackets show performance with modular drying activated, Note: the indicated times are based on an average load observing the indicated loads

Note: times do not take the pre-heat time into account Note: times may vary depending on the load and power supply

## **Accessories**

## A complete array of accessories to expand the Supreme function range.

### **External printer**

Connected to Supreme sterilizers, lets you print cycle-related data on thermal paper or on labels in Barcode or QR code format. To optimise management, the printer can be shared by multiple sterilizers (where enabled) and connected to the same network.

### Front filling kit

This kit allows the sterilizer to be filled frontally via a quick coupling.

### Automatic filling kit

Consists of an external pump powered by the sterilizer, allows demineralised water to be drawn from an external recipient.

#### EV aux kit

This kit allows the sterilizer to be interfaced with surgery demineralising systems so the sterilizer only draws water from the system as and when necessary.


#### Modular tray holder

Modular housing system to load the sterilizer. Supplied with three pairs of shelves to house trays: lets you make full use of the entire chamber diameter and insert even the bulkiest loads.

#### Barcode/QR code reader.

The reader - which can be used on most PCs - lets you complete the traceability process and pair Barcode or QR code labels with the patient.





Ceffa s.c. reserves the right to modify the products illustrated in the photos without prior notice.

According to the standards in force, in extra-EU areas the availability and specifications of some products and/or characteristics may vary. Please contact your local distributor for further information.







**BU Medical Equipment** Sede Legale ed Amministrativa Headquarters

CEFLA s.c.

Via Selice Provinciale, 23/a 40026 Imola (BO) - Italy Tel. +39 0542 653111 Fax +39 0542 653344

#### Stabilimento / Plant

Via Bicocca, 14/C 40026 Imola (BO) - Italy Tel. +39 0542 653441 Fax +39 0542 653601 www.mocom.it

